

# MODIFIED TWIN BLOCK™

Anti-Snoring • Bruxing • Clenching • TMJD



- ✓ Effective
- ✓ Comfortable
- ✓ Affordable
- ✓ Discreet
- ✓ Freedom for Speech
- ✓ Jaw Movement during REM
- ✓ Does not alter the Occlusion
- ✓ Triangular Delta Clasps - optimum retention/minimal fracture

**\$270**  
Shipping  
Included!



*\*Add to the face of the Mx blocks using lightcure Triad for additional Md advancement.*

Serving the dental profession since 1972

Courier service provided by

**Purolator**  
1-888-744-7123

Ship COLLECT under Acct#3575071

## RIDEAU ORTHODONTIC MFG. LTD.

69 Beckwith Street North • Smiths Falls, ON K7A 2B6

1-800-267-7982 • [www.orthodontic.ca](http://www.orthodontic.ca) • [orthodontic@orthodontic.ca](mailto:orthodontic@orthodontic.ca)

### Construction Bite (with patient sitting or standing): *The patient's comfort/balance is an important determining factor.*

- An incisal opening of 2-3mm is required to accommodate the lower anterior capping of acrylic

#### Mandibular Advancement:

**Dental/Skeletal Class I:** Generally 3mm of mandibular advancement along with 6-8mm of vertical opening in the lower 2nd bicuspid region.

**Dental/Skeletal Class II:** Advancement of the mandible to an incisal edge to edge relationship along with 6-8mm of vertical opening in the lower 2nd bicuspid region.

- Comparison of a lateral ceph in CR to the construction bite position may allow you to document the improvement to the airway.
- The construction bite position may open the airway sufficiently while dropping the head of the condyles down/forward in their fossae providing healing space for the TM joint/disc. This position could eliminate/reduce SNORING, bruxing /clenching during sleep.



**CAUTION:** Sleep Apnea can be a life threatening condition. DO NOT attempt to treat a Sleep Apnea patient with any dental oral appliance without instructions by an MD in Sleep Disorders Medicine.

#### Patients with weak musculature:

Hooks are incorporated to allow placement of 5/16" 2 ½ ounce inter-arch elastics. The elastics are generally required for the first three months.

The Modified Twin Block™ should be worn one hour prior to sleep for the first two weeks for patient acceptance.



#### TWIN BLOCK™ FUNCTIONAL THERAPY

Applications in Dentofacial Orthopaedics  
by Dr. William J Clark

Professor Dr. T. M. Graber - Chapter 20  
"Adult Treatment" Foreword

*"It is well known that forward posturing of the mandible has a beneficial effect in most TMD patients."*